

TruTaq DNA Polymerase, 5 U/μL

FB02001020 200 Units

FB02001050 500 Units

Store at -25°C to -15°C

TruTaq DNA Polymerase is an enhanced Taq DNA polymerase optimized for all standard PCR applications, ensuring higher sensitivity, longer PCR products, and higher yields compared to conventional Taq DNA polymerase. TruTaq DNA Polymerase generates PCR products with 3'-dA overhangs. The enzyme is inhibited by dUTP but can incorporate modified nucleotides. Suitable applications include:

- Routine PCR amplification of DNA fragments up to 6 kb from genomic DNA and up to 20 kb from viral DNA.
- RT-PCR.
- Generation of PCR products for TA cloning.

Kit Contents

Reagents	200 Units	500 Units	Description
TruTaq DNA Polymerase, 5 U/μL	40 μL	100 μL	-
10X TruTaq Buffer	1.25 mL	2 x 1.25 mL	Contains KCl and (NH ₄) ₂ SO ₄ at a ratio optimized for robust performance, also includes 20 mM MgCl ₂ .

PCR Reaction Setup:

For multiple reactions, prepare a master mix of components common to all reactions to minimize pipetting error, then dispense appropriate volumes into individual PCR tube before adding template DNA.

Component	50 μL Reaction	Final Concentration
Water, nuclease-free	to 50 μL	-
10X TruTaq Buffer	5 μL	1X
10 mM dNTP Mix	1 μL	0.2 mM each
Forward Primer	varies	0.1-1.0 μM*
Reverse Primer	varies	0.1-1.0 μM*
Template DNA	varies	0.1-1 μg (gDNA) or 0.01-1 ng (plasmid DNA)
TruTaq DNA Polymerase	0.25 μL	1.25 U

* For degenerate primers and primers used for long PCR, higher primer concentrations in the range of 0.3-1 M is recommended.

Thermal Cycling Conditions on Thermal Cycler:

Step	Temp.	Time	Recommendation
Initial denaturation	95°C	1-3 min	Extend to 10 min for GC-rich templates
Denature	95°C	30 sec	Extend to 3-4 min for GC-rich templates
Anneal	Tm-5°C	30 sec	Optimize if non-specific products appear
			Extend 1 minute/kb for longer product > 2 kb
25-40 PCR cycles			Reduce temperature to 68°C for templates >6 kb
Extend	72°C	1 min	
Final extension	72°C	5-15 min	Extend to 30 minutes for TA cloning to ensure complete 3'-dA tailing of the PCR product

Use your PCR product immediately in downstream applications, or store it at -20°C.