

TruFidelis Core PCR Master Mix with HF Buffer

FB02009010 100 Reactions

FB02009050 500 Reactions

Store at -25°C to -15°C

TruFidelis Core DNA Polymerase is a high-fidelity, proofreading DNA polymerase combining a novel enzyme with a processivity-enhancing domain. It offers high performance for all major PCR applications. It generates long amplicons with accuracy and speed, even on challenging templates. The high fidelity makes the TruFidelis Core DNA Polymerase a superior choice for cloning.

- > 50X fidelity compared to Taq polymerase.
- 5'→3' DNA polymerase activity.
- 3'→5' exonuclease activity.
- Generates blunt-end products.
- Amplifies long amplicons such as 7.5 kb genomic and 20 kb λ DNA.

Kit Contents

Reagents	100 Reactions	500 Reactions	Description
2X TruFidelis Core Master Mix with HF Buffer	2 x 1.25 mL	10 x 1.25 mL	Provides 1.5 mM MgCl ₂ at 1X concentration
DMSO	500 μ L	2 x 500 μ L	Recommended for GC-rich amplicons. DMSO is not recommended for amplicons with very low GC% or amplicons that are > 20 kb.

PCR Reaction Setup:

For multiple reactions, prepare a master mix of components common to all reactions to minimize pipetting error, then dispense appropriate volumes into individual PCR tube before adding template DNA.

Component	20 μ L Reaction	50 μ L Reaction	Final Concentration
Water, nuclease-free	to 20 μ L	to 50 μ L	-
2X TruFidelis Core Master Mix with HF Buffer	10 μ L	25 μ L	1X
Forward Primer	varies	varies	0.5 μ M*
Reverse Primer	varies	varies	0.5 μ M*
Template DNA	varies	varies	50-250 ng (gDNA) or 1 pg-10 ng (plasmid DNA) per 50 μ L reaction
DMSO, optional	0.6 μ L	1.5 μ L	3%

*The recommendation for final primer concentration is 0.5 μ M, but it can be varied in a range of 0.2-1.0 μ M, if needed.

Thermal Cycling Conditions on Thermal Cycler:

Step	2-step protocol (for primers with T _m values \geq 69°C)		3-step protocol	
	Temp.	Time	Temp.	Time
Initial denaturation	98°C	30 sec	98°C	30 sec
25-35 PCR cycles	Denature	98°C	5-10 sec	98°C
	Anneal	72°C	15-30 sec/kb	varies
	Extend	72°C	5-10 min	10-30 sec
Final extension	72°C	5-10 min	72°C	15-30 sec/kb
	4°C	hold	4°C	5-10 min
			4°C	hold

Use your PCR product immediately in downstream applications, or store it at -20°C.