

TruTaq Green PCR Master Mix

FB02002020 200 Reactions

Store at -25°C to -15°C

TruTaq Green PCR Master Mix is a ready-to-use solution containing TruTaq DNA polymerase, optimized TruTaq Green buffer, MgCl₂, and dNTPs. It includes a density reagent and two tracking dyes for monitoring electrophoresis progress: The blue dye migrates with 3-5 kb DNA fragments in a 1% agarose gel and the yellow dye migrates faster than 10 bp DNA fragments in 1% agarose gel. The dyes have absorption peaks at 424 nm and 615 nm.

TruTaq Green PCR Master Mix retains all the features of TruTaq DNA Polymerase. It is compatible with downstream workflows such as DNA sequencing, ligation and restriction digestion. Suitable applications include:

- Routine PCR amplification of DNA fragments up to 6 kb from genomic DNA and up to 20 kb from viral DNA.
- RT-PCR.
- Genotyping.
- Generation of PCR products for TA cloning.

Kit Contents

Reagents	200 Reactions	Description
2X TruTaq Green PCR Master Mix	4 x 1.25 mL	Contains dATP, dCTP, dGTP and dTTP, 0.4 mM each, and 4 mM MgCl ₂ .
Water, nuclease-free	4 x 1.25 mL	-

PCR Reaction Setup:

For multiple reactions, prepare a master mix of components common to all reactions to minimize pipetting error, then dispense appropriate volumes into individual PCR tube before adding template DNA.

Component	50 µL Reaction	Final Concentration
Water, nuclease-free	to 50 µL	-
2X TruTaq Green PCR Master Mix	25 µL	1X
Forward Primer	varies	0.1-1.0 µM*
Reverse Primer	varies	0.1-1.0 µM*
Template DNA	varies	0.1-1 µg (gDNA) or 0.01-1 ng (plasmid DNA)

* For degenerate primers and primers used for long PCR, higher primer concentrations in the range of 0.3-1 µM is recommended.

Thermal Cycling Conditions on Thermal Cycler:

Step	Temp.	Time	Recommendation	
Initial denaturation	95°C	1-3 min	Extend to 10 min for GC-rich templates	
25-40 PCR cycles	Denature	95°C	30 sec	Extend to 3-4 min for GC-rich templates
	Anneal	T _m -5°C	30 sec	Optimize if non-specific products appear
	Extend	72°C	1 min	Extend 1 minute/kb for longer product >2 kb Reduce temperature to 68°C for templates >6 kb
Final extension	72°C	5-15 min	Extend to 30 minutes for TA cloning to ensure complete 3'-dA tailing of the PCR product	

Load 5-15 µL of PCR mixture directly on a gel.