

Contents

1. Description
 - 1.1 Principle of the MACS[®] Separation
 - 1.2 Background information
 - 1.3 Applications
 - 1.4 Reagent and instrument requirements
2. Protocol
 - 2.1 Sample preparation
 - 2.2 Magnetic labeling of non-CD4⁺ T cells and fluorescent labeling of CD25⁺ cells
 - 2.3 Magnetic separation: Depletion of non-CD4⁺ T cells
 - 2.4 Magnetic labeling of CD25⁺ cells
 - 2.5 Magnetic separation: Positive selection of CD4⁺CD25⁺ regulatory T cells
3. Example of a separation using the CD4⁺CD25⁺ Regulatory T Cell Isolation Kit
4. References

Warnings

Reagents contain sodium azide. Under acidic conditions sodium azide yields hydrazoic acid, which is extremely toxic. Azide compounds should be diluted with running water before discarding. These precautions are recommended to avoid deposits in plumbing where explosive conditions may develop.

1. Description

This product is for research use only.

Components	<p>1 mL CD4⁺CD25⁺ Regulatory T Cell Biotin-Antibody Cocktail, mouse: Cocktail of biotin-conjugated monoclonal anti-mouse antibodies against CD8a, CD11b, CD45R, CD49b, and Ter-119.</p> <p>2 mL Anti-Biotin MicroBeads: MicroBeads conjugated to monoclonal anti-biotin antibody (isotype: mouse IgG1).</p> <p>0.2 mL CD25-PE, mouse: Monoclonal anti-mouse CD25 antibody conjugated to R-Phycoerythrin (PE) (isotype: rat IgM).</p> <p>1 mL Anti-PE MicroBeads: MicroBeads conjugated to monoclonal anti-PE antibodies (isotype: mouse IgG1).</p>
-------------------	--

Size For 10⁹ total cells, up to 100 separations.

Product format All components are supplied in buffer containing stabilizer and 0.05% sodium azide.

Storage

Store protected from light at 2–8 °C. Do not freeze. The expiration date is indicated on the vial label.

1.1 Principle of the MACS[®] Separation

The isolation of the mouse CD4⁺CD25⁺ regulatory T cells is performed in a two-step procedure. First, the non-CD4⁺ T cells are indirectly magnetically labeled with a cocktail of biotin-conjugated antibodies, as primary labeling reagent, and anti-biotin monoclonal antibodies conjugated to MicroBeads, as secondary labeling reagent. In between the two labeling steps no washing steps are required. In parallel, the cells are labeled with CD25-PE. The magnetically labeled cells are subsequently depleted by separation over a MACS[®] Column, which is placed in the magnetic field of a MACS Separator. The magnetically labeled non-CD4⁺ T cells are retained in the column, while the unlabeled CD4⁺ T cells run through.

In the second step, the CD25⁺ PE-labeled cells are magnetically labeled with Anti-PE MicroBeads and isolated by positive selection from the pre-enriched CD4⁺ cell fraction by separation over a MACS Column, which is placed in the magnetic field of a MACS Separator. The magnetically labeled CD4⁺CD25⁺ cells are retained in the column, while the unlabeled CD4⁺CD25⁻ cells run through and can be used for suppression assays as responder T cells.

After removing the column from the magnetic field, the magnetically retained CD4⁺CD25⁺ cells can be eluted as the positively selected cell fraction. To increase the purity, the positively selected cell fraction containing the CD4⁺CD25⁺ cells must be separated over a second column.

Mouse splenocytes: Depletion of non-CD4⁺ T cells

1. Indirect magnetic labeling of non-CD4⁺ T cells with CD4⁺CD25⁺ Regulatory T Cell Biotin-Antibody Cocktail and Anti-Biotin MicroBeads.
2. Labeling of CD25⁺ cells with CD25-PE.
3. Magnetic separation using an LD Column or an autoMACS Column (program "Depl025").

CD4⁺ T cells (flow-through fraction): Positive selection of CD4⁺CD25⁺ cells

1. Indirect magnetic labeling of CD25⁺ cells with Anti-PE MicroBeads.
2. Magnetic separation using two MS Columns or an autoMACS Column (program "Posseld2").

CD4⁺CD25⁺ regulatory T cells (eluted fraction) CD4⁺CD25⁻ T cells (flow-through fraction)

1.2 Background information

The CD4⁺CD25⁺ Regulatory T Cell Isolation Kit has been developed for the isolation of mouse CD4⁺CD25⁺ regulatory T cells from single-cell suspensions of spleen or lymph nodes without centrifugational washing steps during the separation.

CD4⁺CD25⁺ immunoregulatory T cells have been shown to actively suppress immune responses against autologous and foreign antigens *in vivo* and *in vitro*. CD25, the IL-2R α chain, is also expressed on activated CD8⁺ T cells, dendritic cells (DCs), and B cells.

The kit contains a cocktail of lineage specific biotin-conjugated antibodies against CD8, CD11b, CD45R, CD49b, Ter-119, and Anti-Biotin MicroBeads for depletion of non-CD4⁺ T cells, as well as CD25-PE and Anti-PE MicroBeads for subsequent positive selection of CD4⁺CD25⁺ regulatory T cells.

1.3 Applications

- Isolation of CD4⁺CD25⁺ regulatory T cells from single-cell suspensions of spleen or lymph nodes for:
 - Treg expansion by using, e.g., the Treg Expansion Kit, mouse (# 130-095-925),
 - *in vitro* Treg suppression assays in combination with CD4⁺CD25⁻ cell fraction as T responder/effector cells,
 - co-culture experiments with DCs to study priming of DCs for tolerance induction *in vitro* and after adoptive transfer of primed DCs *in vivo*,¹
 - adoptive transfer experiments (e.g. from UV-exposed mice) to analyze the role of regulatory T cells during induction and elicitation of hapten-specific tolerance.²

1.4 Reagent and instrument requirements

- Buffer: Prepare a solution containing phosphate-buffered saline (PBS), pH 7.2, 0.5% bovine serum albumin (BSA), and 2 mM EDTA by diluting MACS BSA Stock Solution (# 130-091-376) 1:20 with autoMACS[®] Rinsing Solution (# 130-091-222). Keep buffer cold (2–8 °C). Degas buffer before use, as air bubbles could block the column.
 - ▲ **Note:** EDTA can be replaced by other supplements such as anticoagulant citrate dextrose formula-A (ACD-A) or citrate phosphate dextrose (CPD). BSA can be replaced by other proteins such as mouse serum albumin, mouse serum, or fetal bovine serum (FBS). Buffers or media containing Ca²⁺ or Mg²⁺ are not recommended for use.
- MACS Columns and MACS Separators: Depletion of non-CD4 T cells can be performed on an LD Column. The subsequent positive selection of CD4⁺CD25⁺ cells can be performed on two MS Columns. Positive selection and depletion can also be performed by using the autoMACS Pro Separator.

Column	Max. number of labeled cells	Max. number of total cells	Separator
Positive selection			
MS	10 ⁷	2×10 ⁸	MiniMACS, OctoMACS, SuperMACS II
Depletion			
LD	10 ⁸	5×10 ⁸	MidiMACS, QuadroMACS, SuperMACS II
Positive selection or depletion			
autoMACS	2×10 ⁸	4×10 ⁹	autoMACS Pro

▲ **Note:** Column adapters are required to insert certain columns into the SuperMACS[™] II Separators. For details refer to the respective MACS Separator data sheet.

- (Optional) Fluorochrome-conjugated antibodies for flow cytometric analysis, e.g., CD4 Antibody, anti-mouse, FITC. For more information about antibodies refer to www.miltenyibiotec.com/antibodies.
- (Optional) FoxP3 Antibody, anti-mouse, APC and FoxP3 Staining Buffer Set (# 130-093-142).
- (Optional) Propidium Iodide Solution (# 130-093-233) or 7-AAD Staining Solution (# 130-111-568) for flow cytometric exclusion of dead cells.
- (Optional) Dead Cell Removal Kit (# 130-090-101) for the depletion of dead cells.
- (Optional) Pre-Separation Filters (30 μm) (# 130-041-407) to remove cell clumps.

2. Protocol

2.1 Sample preparation

Prepare a single-cell suspension from spleen or lymph nodes using manual methods or the gentleMACS[™] Dissociators.

For details refer to www.gentleMACS.com/protocols.

▲ **Note:** The kit is not optimized for the isolation of regulatory T cells from blood and thymus.

▲ **Note:** Red blood cell lysis or density gradient centrifugation is not necessary, since the CD4⁺CD25⁺ T Cell Biotin-Antibody Cocktail contains Anti-Ter-119 antibody.



2.2 Magnetic labeling of non-CD4⁺ cells and fluorescent labeling of CD25⁺ cells

▲ Work fast, keep cells cold, and use pre-cooled solutions. This will prevent capping of antibodies on the cell surface and non-specific cell labeling.

▲ Volumes for magnetic labeling given below are for up to 10⁷ total cells. When working with fewer than 10⁷ cells, use the same volumes as indicated. When working with higher cell numbers, scale up all reagent volumes and total volumes accordingly (e.g. for 2×10⁷ total cells, use twice the volume of all indicated reagent volumes and total volumes).

▲ For optimal performance it is important to obtain a single-cell suspension before magnetic labeling. Pass cells through 30 μm nylon mesh (Pre-Separation Filters (30 μm), # 130-041-407) to remove cell clumps which may clog the column. Moisten filter with buffer before use.

▲ The recommended incubation temperature is 2–8 °C. Higher temperatures and/or longer incubation times may lead to non-specific cell labeling. Working on ice may require increased incubation times.

1. Determine cell number.
2. Centrifuge cell suspension at 300×g for 10 minutes. Aspirate supernatant completely.
3. Resuspend cell pellet in 40 μL of buffer per 10⁷ total cells.
4. Add 10 μL of the CD4⁺CD25⁺ Regulatory T Cell Biotin-Antibody Cocktail per 10⁷ total cells.
5. Mix well and incubate for 10 minutes in the refrigerator (2–8 °C).

- Add 38 μL of buffer, 20 μL of Anti-Biotin MicroBeads, and 2 μL of CD25-PE antibody per 10^7 total cells.

▲ **Note:** CD4⁺CD25⁺ target cells are always stained with CD25-PE due to the separation procedure and can be detected via flow cytometry.

- Mix well and incubate for additional 15 minutes in the refrigerator (2–8 °C).

- Proceed to magnetic separation (2.3).

▲ **Note:** A minimum of 500 μL is required for magnetic separation. If necessary, add buffer to the cell suspension.



2.3 Magnetic separation: Depletion of non-CD4⁺ T cells

▲ Always wait until the column reservoir is empty before proceeding to the next step.

Depletion with LD Columns

- Place LD Column in the magnetic field of a suitable MACS Separator. For details refer to the LD Column data sheet.
- Prepare column by rinsing with 2 mL of buffer.
- Apply cell suspension onto the column.
- Collect unlabeled cells that pass through and wash column with 2 \times 1 mL of buffer. Collect total flow-through; this is the unlabeled CD4⁺ cell fraction. Perform washing steps by adding buffer two times. Only add new buffer when the column reservoir is empty.
- Proceed to 2.4 for the magnetic labeling of CD25⁺ T cells.

Depletion with the autoMACS Pro Separator

▲ Refer to the respective user manual for instructions on how to use the autoMACS Pro Separator.

▲ Buffers used for operating the autoMACS Pro Separator should have a temperature of ≥ 10 °C.

- Prepare and prime the instrument.
- Apply tube containing the sample and provide tubes for collecting the labeled and unlabeled cell fractions. Place sample tube in row A of the tube rack and the fraction collection tubes in rows B and C.
- For a standard separation choose the following program:

Depletion: Depl025

Collect negative fraction in row B of the tube rack.

- Proceed to 2.4 for the magnetic labeling of CD25⁺ cells.



2.4 Magnetic labeling of CD25⁺ cells

▲ Volumes for magnetic labeling given below are for an initial starting cell number of up to 10^7 total cells. For higher initial cell numbers, scale up all volumes accordingly.

- Centrifuge cell suspension at 300 \times g for 5 minutes. Aspirate supernatant completely.
- Resuspend cell pellet in 90 μL of buffer.
- Add 10 μL of Anti-PE MicroBeads.
- Mix well and incubate for 15 minutes in the refrigerator (2–8 °C).
- Proceed to magnetic separation (2.5).

▲ **Note:** A minimum of 500 μL is required for magnetic separation. If necessary, add buffer to the cell suspension.



2.5 Magnetic separation: Positive selection of CD4⁺CD25⁺ regulatory T cells

Positive selection with MS Columns

▲ To achieve highest purities, perform two consecutive column runs.

- Place MS Column in the magnetic field of a suitable MACS Separator. For details refer to the MS Column data sheet.
 - Prepare column by rinsing with 500 μL of buffer.
 - Apply cell suspension onto the column. Collect flow-through containing CD4⁺CD25⁻ cells.
 - Wash column with 2 \times 500 μL of buffer. Collect unlabeled cells that pass through and combine with the flow-through from step 3.
- ▲ **Note:** Perform washing steps by adding buffer aliquots as soon as the column reservoir is empty.
- Remove column from the separator and place it on a suitable collection tube.
 - Pipette 1 mL of buffer onto the column. Immediately flush out the magnetically labeled cells (CD4⁺CD25⁺ cells) by firmly pushing the plunger into the column.
- ▲ **Note:** To perform a second column run, you may elute the cells directly from the first onto the second, equilibrated column instead of a collection tube.
- Repeat the magnetic separation procedure as described in steps 1 to 6 by using a new MS Column.

Positive selection with the autoMACS Pro Separator


- Prepare and prime the instrument.
- Apply tube containing the sample and provide tubes for collecting the labeled and unlabeled cell fractions. Place sample tube in row A of the tube rack and the fraction collection tubes in rows B and C.
- For a standard separation choose the following program:

Positive selection: Posseld2

Collect positive fraction in row C of the tube rack. This is the enriched CD4⁺CD25⁺ T cell fraction.

3. Example of a separation using the CD4⁺CD25⁺ Regulatory T Cell Isolation Kit

CD4⁺CD25⁺ regulatory T cells were isolated from mouse spleen cell suspension by using the CD4⁺CD25⁺ Regulatory T Cell Isolation Kit, an LD and two MS Columns, a MidiMACS™ Separator and a MiniMACS™ Separator. The cells were fluorescently stained with CD25-PE and CD4-FITC (A) or Anti-FoxP3-APC (B) and analyzed by flow cytometry using the MACSQuant® Analyzer. Cell debris and dead cells were excluded from analysis based on scatter signals and propidium iodide fluorescence.



4. References

1. Fallarino, F. *et al.* (2003) Modulation of tryptophan catabolism by regulatory T cells. *Nat. Immunol.* 4(12): 1206–1212.
2. Schwarz, A. *et al.* (2004) Ultraviolet radiation-induced regulatory T cells not only inhibit the induction but can suppress the effector phase of contact hypersensitivity. *J. Immunol.* 172: 1036–1043.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

autoMACS, gentleMACS, MACS, MACSQuant, MidiMACS, the Miltenyi Biotec logo, MiniMACS, OctoMACS, QuadroMACS, and SuperMACS are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2022 Miltenyi Biotec and/or its affiliates. All rights reserved.