

## 10X TBE Electrophoresis Buffer

FB02511100 1 L

Stored at room temperature

### Kit Contents

Reagents	Volume	1X Buffer Composition	Notes
10X TBE Electrophoresis Buffer	1 L	89 mM Tris, 89 mM boric acid, 2 mM EDTA	<ul style="list-style-type: none"><li>• Dilute to working concentration of 1X before use</li><li>• Use fresh 1X buffer for each electrophoresis</li></ul>

### Applications

- Nucleic acid agarose and polyacrylamide gel electrophoresis.
- Agarose and polyacrylamide gel preparation.