

10mM dNTP Mix

Molecular biology grade

FB02512001 1 mL

FB02512005 5 x 1mL

Store at -25°C to -15°C

Description

dNTP Mix contains aqueous solution of dATP, dCTP, dGTP and dTTP, each at a final concentration of 10 mM, which offers the possibility to reduce the number of pipetting steps, and the risk of reaction set up errors.

Applications

For use in PCR, real-time PCR, high fidelity and long PCR, LAMP-PCR, cDNA synthesis, RT-PCR, RDA, MDA, DNA labeling, and DNA sequencing.